Anleitung

Dr. Arabin GmbH & Co KG
Alfred-Herrhausen-Str. 44, 58455 Witten
Tel: +49 2302 189214 Fax: +49 2302 189216
www.dr-arabin.de info@dr-arabin.de

Instructions

Dr. Arabin GmbH & Co KG
Alfred-Herrhausen-Str. 44, 58455 Witten
Tel: (49) 0 2302 189214 Fax: (49) 0 2302 189216
www.dr-arabin.de info@dr-arabin.de

Cerclage Pessar

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D-CA+ASQ

(Type A without perforations / Type ASQ perforated)

Application: The cervical pessary was designed, to prevent spontaneous preterm birth for pregnant women who are exposed to physical stress (long periods of standing or walking), increased intrauterine pressure (polyhydramnios, multiple pregnancies) and/ or early clinical and/or ultrasound signs of a short cervix.


The pessary consists of biocompatible silicone, it can be folded and then inserted without pain. Perforated pessaries (ASQ) allow a better passage of discharge, which is increased during treatment, but usually not combined with a change of the bacterial pattern or infection. Nevertheless, it may be indicated to perform a vaginal swab to know the kind of vaginal colonization in patients at increased risk for preterm birth.

The pessary placement should start at an early stage of precocious cervical ripening. The pessary is not meant to close the cervix completely but rather to reduce cervical stretching and to change the utero-cervical angle.

Sizes: Cervical pessaries are indicated according to

- a) the outer diameter (65 mm for nulliparous or small women and 70mm for multiparous or tall women),
- b) the height of the curvature (17mm rarely during the first trimester, 21 mm for smaller patients and singleton pregnancy, 25mm for taller patients and multiple pregnancy, 30mm: rarely in women with genital prolapse),
- c) the inner diameter (32 mm for average cervical shortening and 35 mm for an unusually wide cervix or U-shaped funneling to avoid further manipulation) (www.dr-arabin.de). The recommendations are relative and can be individually modified.

Use: The obstetrician places the cervical pessary with the patient is in a lithotomy position. It is vital that **the curvature of the pessary with the smaller diameter shows upwards**. During insertion, the pessary stays folded until the upper vaginal fornix is reached and is then gently "pushed" as high as possible into the posterior vaginal vault with the smaller diameter around the complete cervix. The anterior rim may be gently pushed towards the sacrum (video, www.dr-arabin.de). Some patients indicate "relief of prolapse complaints" after placement when they get up. In any case, the patient should not any longer feel the pessary.

In asymptomatic patients, the pessary should be removed before delivery at around 37 weeks or before a planned Cesarean. In rare cases (discomfort, intermediate blood loss) the pessary may be removed and another size adapted. In general, manipulations are better avoided. In patients with premature rupture of membranes and signs of chorioamnionitis, persistent vaginal blood loss and regular or painful contractions a speculum examination should be performed and the pessary should be removed to avoid cervical lesions or ascending infection. If the removal seems difficult (e.g. in rare cases with prolapse or contractions) the pessary should be cut with atraumatic scissors instead of using violence.

Since the visualization of the cervix by transvaginal sonography may be impaired by the pessary itself, the transvaginal transducer should be directed "on top" of the upper cervical rim placing it onto the anterior cervix to avoid manipulation. Alternatively it can be visualized from the posterior vaginal fornix. Usually, the position of the (rotated) pessary can easily be controlled manually without touching the cervix.

Side Effects/Complications: The placement of the pessary with the smaller diameter upwards and some angulation is essential. Vaginal discharge is a frequent side effect. However, the composition of the vaginal flora is usually not be substantially altered and no antibiotics are needed only to treat discharge. It is essential not to confuse discharge with premature rupture of membranes. In case of uncertainty, one should perform transabdominal ultrasound and biochemical tests before removing the pessary unnecessarily.

In some women – mainly if they might have some degree of genital prolapse- the pessary might descend during standing or pushing towards the introitus (with the cervix) and patients or physicians might confuse this with protruding membranes. Then, patients need to be re-assured and the pessary should be pushed upwards provided that premature rupture of membranes is excluded.

This is a **PRODUCT FOR SINGLE USE ONLY**. In case of repetitive use an infection can not be excluded. It may be stored by room temperature between 1 and 30 °C in the absence of UV irradiation. In addition, the direct contact to reactive and non-polar media, e.g. gases like ozone or solvents like mineral oils should be avoided.

GB-CA+ASQ

Instructions

Dr. Arabin GmbH & Co KG
 Alfred-Herrhausen-Str. 44, 58455 Witten
 Tel: (49) 0 2302 189214 Fax: (49) 0 2302 189216
www.dr-arabin.de info@dr-arabin.de

Pessaire-cerclage

(Type A sans perforation / Type ASQ perforé)

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Dr. Arabin GmbH & Co KG im FEZ